

Product Name: GABAA R α 4 Rabbit Polyclonal Antibody
Catalog #: APRab11234



Summary

Production Name	GABAA R α 4 Rabbit Polyclonal Antibody
Description	Rabbit Polyclonal Antibody
Host	Rabbit
Application	WB,ELISA
Reactivity	Human,Mouse,Rat

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Polyclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% New type preservative N.
Purification	Affinity purification

Immunogen

Gene Name	GABRA4
Alternative Names	GABRA4; Gamma-aminobutyric acid receptor subunit alpha-4; GABA(A) receptor subunit alpha-4
Gene ID	2557.0
SwissProt ID	P48169.The antiserum was produced against synthesized peptide derived from human GABRA4. AA range:81-130

Application

Dilution Ratio	WB 1:500 - 1:2000. ELISA: 1:40000. Not yet tested in other applications.
Molecular Weight	60kD

Product Name: GABAA R α 4 Rabbit Polyclonal Antibody
Catalog #: APRab11234




Background

Gamma-aminobutyric acid (GABA) is the major inhibitory neurotransmitter in the mammalian brain where it acts at GABA-A receptors, which are ligand-gated chloride channels. Chloride conductance of these channels can be modulated by agents such as benzodiazepines that bind to the GABA-A receptor. At least 16 distinct subunits of GABA-A receptors have been identified. This gene encodes subunit alpha-4, which is involved in the etiology of autism and eventually increases autism risk through interaction with another subunit, gamma-aminobutyric acid receptor beta-1 (GABRB1). Alternatively spliced transcript variants encoding different isoforms have been found in this gene.[provided by RefSeq, Feb 2011],function:GABA, the major inhibitory neurotransmitter in the vertebrate brain, mediates neuronal inhibition by binding to the GABA/benzodiazepine receptor and opening an integral chloride channel.,induction:The alpha4 beta2 gamma 2L receptor is not repressed by diazepam.,online information:Forbidden fruit - Issue 56 of March 2005,similarity:Belongs to the ligand-gated ionic channel (TC 1.A.9) family.,subunit:Generally pentameric. There are five types of GABA(A) receptor chains: alpha, beta, gamma, delta, and rho.,


Research Area

Neuroactive ligand-receptor interaction;

Image Data




Western blot analysis of GABRA4 Antibody. The lane on the right is blocked with the GABRA4 peptide.



Western blot analysis of the lysates from COLO205 cells using GABRA4 antibody.

Product Name: GABAA R α 4 Rabbit Polyclonal Antibody
Catalog #: APRab11234



Western Blot analysis of various cells using GABAA R α 4 Polyclonal Antibody diluted at 1: 1000

Note

For research use only.