

Product Name: Galactosidase β 1L3 Rabbit Polyclonal Antibody
Catalog #: APRab11262

Summary

Production Name	Galactosidase β 1L3 Rabbit Polyclonal Antibody
Description	Rabbit Polyclonal Antibody
Host	Rabbit
Application	IHC,ELISA
Reactivity	Human,Rat,Mouse

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Polyclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% New type preservative N.
Purification	Affinity purification

Immunogen

Gene Name	GLB1L3
Alternative Names	GLB1L3; Beta-galactosidase-1-like protein 3
Gene ID	112937.0
SwissProt ID	Q8NCI6.The antiserum was produced against synthesized peptide derived from human BGLL3. AA range:291-340

Application

Dilution Ratio	IHC 1:100-1:300 ELISA: 1:40000
Molecular Weight	

Background


Product Name: Galactosidase β 1L3 Rabbit Polyclonal Antibody
Catalog #: APRab11262



similarity:Belongs to the glycosyl hydrolase 35 family.,similarity:Belongs to the glycosyl hydrolase 35 family.,

Research Area

Image Data



Immunohistochemistry analysis of paraffin-embedded human colon carcinoma, using GLB1L3 Antibody. The picture on the right is blocked with the synthesized peptide.

Note

For research use only.