

Product Name: IgD Rabbit Polyclonal Antibody
Catalog #: APRab12419



Summary

Production Name	IgD Rabbit Polyclonal Antibody
Description	Rabbit Polyclonal Antibody
Host	Rabbit
Application	IHC,ELISA
Reactivity	Human,Rat,Mouse

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Polyclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% New type preservative N.
Purification	Affinity purification

Immunogen

Gene Name	IGHD
Alternative Names	Ig delta chain C region
Gene ID	
SwissProt ID	P01880.Synthetic peptide from human protein at AA range: 111-160

Application

Dilution Ratio	IHC 1:50-200 ELISA 1:10000-20000
Molecular Weight	

Background


similarity:Contains 3 Ig-like (immunoglobulin-like) domains.,similarity:Contains 3 Ig-like (immunoglobulin-like) domains.,

Product Name: IgD Rabbit Polyclonal Antibody
Catalog #: APRab12419



Research Area

Image Data



Immunohistochemical analysis of paraffin-embedded human-tonsil, antibody was diluted at 1:200

Note

For research use only.