

Product Name: IQGAP3 Rabbit Polyclonal Antibody
Catalog #: APRab12723



Summary

Production Name	IQGAP3 Rabbit Polyclonal Antibody
Description	Rabbit Polyclonal Antibody
Host	Rabbit
Application	WB,ELISA
Reactivity	Human,Mouse,Rat

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Polyclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% New type preservative N.
Purification	Affinity purification

Immunogen

Gene Name	IQGAP3
Alternative Names	IQGAP3; Ras GTPase-activating-like protein IQGAP3
Gene ID	128239.0
SwissProt ID	Q86VI3.The antiserum was produced against synthesized peptide derived from human IQGAP3. AA range:107-156

Application

Dilution Ratio	WB 1:500-2000 ELISA 2000-20000
Molecular Weight	180kD

Background

Product Name: IQGAP3 Rabbit Polyclonal Antibody
Catalog #: APRab12723




similarity:Contains 1 CH (calponin-homology) domain.,similarity:Contains 1 Ras-GAP domain.,similarity:Contains 4 IQ domains.,similarity:Contains 1 CH (calponin-homology) domain.,similarity:Contains 1 Ras-GAP domain.,similarity:Contains 4 IQ domains.,


Research Area

Regulates Actin and Cytoskeleton;

Image Data



Western blot analysis of the lysates from 293 cells using IQGAP3 antibody.



Western Blot analysis of various cells using IQGAP3 Polyclonal Antibody diluted at 1: 1000

Note

For research use only.