

**Product Name: MBL-C Rabbit Polyclonal Antibody**  
**Catalog #: APRab13686**



## Summary

<b>Production Name</b>	MBL-C Rabbit Polyclonal Antibody
<b>Description</b>	Rabbit Polyclonal Antibody
<b>Host</b>	Rabbit
<b>Application</b>	WB,IHC,
<b>Reactivity</b>	Human,Mouse

## Performance

<b>Conjugation</b>	Unconjugated
<b>Modification</b>	Unmodified
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>Form</b>	Liquid
<b>Storage</b>	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
<b>Buffer</b>	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% New type preservative N.
<b>Purification</b>	Affinity purification

## Immunogen

<b>Gene Name</b>	MBL2 COLEC1 MBL
<b>Alternative Names</b>	MBL2 COLEC1 MBL
<b>Gene ID</b>	4153.0
<b>SwissProt ID</b>	P11226.Synthetic peptide from human protein at AA range: 21-70

## Application

<b>Dilution Ratio</b>	WB 1:500-2000;IHC 1:50-300
<b>Molecular Weight</b>	27kD

## Background

This gene encodes the soluble mannose-binding lectin or mannose-binding protein found in serum. The protein encoded

**Product Name: MBL-C Rabbit Polyclonal Antibody**  
**Catalog #: AP Rab13686**




belongs to the collectin family and is an important element in the innate immune system. The protein recognizes mannose and N-acetylglucosamine on many microorganisms, and is capable of activating the classical complement pathway. Deficiencies of this gene have been associated with susceptibility to autoimmune and infectious diseases. [provided by RefSeq, Jul 2008],disease:Genetic variations in MBL2 are associated with susceptibility to hepatitis B virus infection (HBV infection) [MIM:610424]. Approximately one third of all cases of cirrhosis and half of all cases of hepatocellular carcinoma can be attributed to chronic HBV infection. HBV infection may result in subclinical or asymptomatic infection, acute self-limited hepatitis, or fulminant hepatitis requiring liver transplantation.,disease:There is an association between low levels of MBL2 and a defect of opsonization which results in susceptibility to frequent and chronic infections.,function: Binds mannose and N-acetylglucosamine in a calcium-dependent manner. Is capable of host defense against pathogens, by activating the classical complement pathway independently of the antibody.,online information:Mannose-binding protein,similarity:Contains 1 C-type lectin domain.,similarity:Contains 1 collagen-like domain.,subunit:Oligomeric complex of 6 set of homotrimers. Interacts with MASP1 and MASP2; the interaction is calcium-dependent.,


## Research Area

Complement and coagulation cascades;

## Image Data




Western Blot analysis of mouse-kidney cells using Antibody diluted at 500. Secondary antibody was diluted at 1:20000




Immunohistochemical analysis of paraffin-embedded human-liver, antibody was diluted at 1:200

**Product Name: MBL-C Rabbit Polyclonal Antibody**  
**Catalog #: APRab13686**


---



Immunohistochemical analysis of paraffin-embedded human-liver, antibody was diluted at 1:200



Immunohistochemical analysis of paraffin-embedded human-brain, antibody was diluted at 1:200



Immunohistochemical analysis of paraffin-embedded human-brain, antibody was diluted at 1:200

## **Note**

For research use only.