

Product Name: MEGF6 Rabbit Polyclonal Antibody
Catalog #: APRab13788



Summary

Production Name	MEGF6 Rabbit Polyclonal Antibody
Description	Rabbit Polyclonal Antibody
Host	Rabbit
Application	WB
Reactivity	Human,Rat,Mouse

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Polyclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Liquid in PBS containing 50% glycerol, and 0.02% New type preservative N.
Purification	Affinity purification

Immunogen

Gene Name	MEGF6 EGFL3 KIAA0815
Alternative Names	
Gene ID	1953.0
SwissProt ID	O75095.Synthesized peptide derived from human protein . at AA range: 40-120

Application

Dilution Ratio	WB 1:500-2000 ELISA 1:5000-20000
Molecular Weight	169kD

Background


similarity:Contains 25 EGF-like domains.,similarity:Contains 25 EGF-like domains.,

Product Name: MEGF6 Rabbit Polyclonal Antibody
Catalog #: APRab13788



Research Area

Image Data



Western Blot analysis of HEK293 lysis, using primary antibody at 1:1000 dilution. Secondary antibody was diluted at 1:10000

Note

For research use only.