

**Product Name: MRP-L12 Rabbit Polyclonal Antibody**  
**Catalog #: APRab14110**



## Summary

<b>Production Name</b>	MRP-L12 Rabbit Polyclonal Antibody
<b>Description</b>	Rabbit Polyclonal Antibody
<b>Host</b>	Rabbit
<b>Application</b>	WB,ELISA
<b>Reactivity</b>	Human,Mouse

## Performance

<b>Conjugation</b>	Unconjugated
<b>Modification</b>	Unmodified
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>Form</b>	Liquid
<b>Storage</b>	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
<b>Buffer</b>	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% New type preservative N.
<b>Purification</b>	Affinity purification

## Immunogen

<b>Gene Name</b>	MRPL12
<b>Alternative Names</b>	MRPL12; RPML12; 39S ribosomal protein L12; mitochondrial; L12mt; MRP-L12; 5c5-2
<b>Gene ID</b>	6182.0
<b>SwissProt ID</b>	P52815.The antiserum was produced against synthesized peptide derived from human MRPL12. AA range:31-80

## Application

<b>Dilution Ratio</b>	WB 1:500-2000;ELISA 2000-20000
<b>Molecular Weight</b>	21kD

## Background


**Product Name: MRP-L12 Rabbit Polyclonal Antibody**  
**Catalog #: APRab14110**



Mammalian mitochondrial ribosomal proteins are encoded by nuclear genes and help in protein synthesis within the mitochondrion. Mitochondrial ribosomes (mitoribosomes) consist of a small 28S subunit and a large 39S subunit. They have an estimated 75% protein to rRNA composition compared to prokaryotic ribosomes, where this ratio is reversed. Another difference between mammalian mitoribosomes and prokaryotic ribosomes is that the latter contain a 5S rRNA. Among different species, the proteins comprising the mitoribosome differ greatly in sequence, and sometimes in biochemical properties, which prevents easy recognition by sequence homology. This gene encodes a 39S subunit protein which forms homodimers. In prokaryotic ribosomes, two L7/L12 dimers and one L10 protein form the L8 protein complex. [provided by RefSeq, Jul 2008],similarity:Belongs to the ribosomal protein L12P family.,

## Research Area

## Image Data



Western blot analysis of lysates from COLO cells, using MRPL12 Antibody. The lane on the right is blocked with the synthesized peptide.



Western Blot analysis of various cells using MRP-L12 Polyclonal Antibody

## Note

For research use only.

**Product Name: MRP-L12 Rabbit Polyclonal Antibody**  
**Catalog #: APRab14110**

---

