

Product Name: NOL4 Rabbit Polyclonal Antibody
Catalog #: APRab14789



Summary

Production Name	NOL4 Rabbit Polyclonal Antibody
Description	Rabbit Polyclonal Antibody
Host	Rabbit
Application	WB
Reactivity	Human,Mouse

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Polyclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Liquid in PBS containing 50% glycerol, and 0.02% New type preservative N.
Purification	Affinity purification

Immunogen

Gene Name	NOL4 NOLP HRIHFB2255
Alternative Names	
Gene ID	8715.0
SwissProt ID	O94818.Synthesized peptide derived from human protein . at AA range: 50-130

Application

Dilution Ratio	WB 1:500-2000 ELISA 1:5000-20000
Molecular Weight	70kD

Background


tissue specificity:Expressed predominantly in fetal brain, adult brain and testis.,tissue specificity:Expressed predominantly in fetal brain, adult brain and testis.,

Product Name: NOL4 Rabbit Polyclonal Antibody
Catalog #: APRab14789



Research Area

Image Data



Western Blot analysis of HEK293 lysis, using primary antibody at 1:1000 dilution. Secondary antibody was diluted at 1:10000

Note

For research use only.