

Product Name: Olfactory receptor 8J1 Rabbit Polyclonal Antibody
Catalog #: APRab15331

Summary

Production Name	Olfactory receptor 8J1 Rabbit Polyclonal Antibody
Description	Rabbit Polyclonal Antibody
Host	Rabbit
Application	IF, WB,
Reactivity	Human, Rat, Mouse

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Polyclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% New type preservative N.
Purification	Affinity purification

Immunogen

Gene Name	OR8J1
Alternative Names	OR8J1; Olfactory receptor 8J1; Olfactory receptor OR11-183
Gene ID	219477.0
SwissProt ID	Q8NGP2. The antiserum was produced against synthesized peptide derived from human OR8J1. AA range: 233-282

Application

Dilution Ratio	WB 1:500 - 1:2000. IF 1:200 - 1:1000. ELISA: 1:10000. Not yet tested in other applications.
Molecular Weight	35kD

Product Name: Olfactory receptor 8J1 Rabbit Polyclonal Antibody
Catalog #: APRab15331




Background

olfactory receptor family 8 subfamily J member 1(OR8J1) Homo sapiens Olfactory receptors interact with odorant molecules in the nose, to initiate a neuronal response that triggers the perception of a smell. The olfactory receptor proteins are members of a large family of G-protein-coupled receptors (GPCR) arising from single coding-exon genes. Olfactory receptors share a 7-transmembrane domain structure with many neurotransmitter and hormone receptors and are responsible for the recognition and G protein-mediated transduction of odorant signals. The olfactory receptor gene family is the largest in the genome. The nomenclature assigned to the olfactory receptor genes and proteins for this organism is independent of other organisms. [provided by RefSeq, Jul 2008],function:Odorant receptor ,similarity:Belongs to the G-protein coupled receptor 1 family.,


Research Area

Olfactory transduction;

Image Data




Immunofluorescence analysis of MCF7 cells, using OR8J1 Antibody. The picture on the right is blocked with the synthesized




Western blot analysis of lysates from Jurkat cells, using OR8J1 Antibody. The lane on the right is blocked with the

Product Name: Olfactory receptor 8J1 Rabbit Polyclonal Antibody
Catalog #: APRab15331



Western blot analysis of the lysates from HepG2 cells using OR8J1 antibody.



Western Blot analysis of Jurkat cells using Olfactory receptor 8J1 Polyclonal Antibody

Note

For research use only.