

Product Name: Pan myristoylation Rabbit Polyclonal Antibody
Catalog #: APRab15722



Summary

Production Name	Pan myristoylation Rabbit Polyclonal Antibody
Description	Rabbit Polyclonal Antibody
Host	Rabbit
Application	WB
Reactivity	Species independent

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Polyclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% New type preservative N.
Purification	Affinity purification

Immunogen

Gene Name	
Alternative Names	
Gene ID	
SwissProt ID	.Synthetic myristoylated peptide

Application

Dilution Ratio	WB 1:500-2000, ELISA 1:10000-20000
Molecular Weight	130 80 35kD


Background

**Product Name: Pan myristoylation Rabbit Polyclonal
Antibody**
Catalog #: APRab15722



Research Area

Image Data



Western Blot analysis of KB cells using Antibody diluted at 500. Secondary antibody was diluted at 1:20000

Note

For research use only.