

**Product Name: Peroxin 11γ Rabbit Polyclonal Antibody**  
**Catalog #: APRab15983**



## Summary

<b>Production Name</b>	Peroxin 11γ Rabbit Polyclonal Antibody
<b>Description</b>	Rabbit Polyclonal Antibody
<b>Host</b>	Rabbit
<b>Application</b>	IHC,ELISA
<b>Reactivity</b>	Human,Rat,Mouse

## Performance

<b>Conjugation</b>	Unconjugated
<b>Modification</b>	Unmodified
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>Form</b>	Liquid
<b>Storage</b>	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
<b>Buffer</b>	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% New type preservative N.
<b>Purification</b>	Affinity purification

## Immunogen

<b>Gene Name</b>	PEX11G
<b>Alternative Names</b>	PEX11G; PEX11C; Peroxisomal membrane protein 11C; Peroxin-11C; Peroxisomal biogenesis factor 11C; Protein PEX11 homolog gamma; PEX11-gamma
<b>Gene ID</b>	92960.0
<b>SwissProt ID</b>	Q96HA9.The antiserum was produced against synthesized peptide derived from human PEX11G. AA range:11-60

## Application

<b>Dilution Ratio</b>	IHC 1:100-1:300 ELISA: 1:20000
<b>Molecular Weight</b>	

## Background


**Product Name: Peroxin 11γ Rabbit Polyclonal Antibody**  
**Catalog #: AP Rab15983**



The protein encoded by this gene is a member of the PEX11 family. This family is reported to regulate the number and size of peroxisomes in evolutionarily distant organisms. The protein encoded by this gene may induce clustering of peroxisomes. Alternative splicing results in multiple transcript variants that encode different protein isoforms. [provided by RefSeq, Jul 2012], similarity: Belongs to the peroxin-11 family.,

## Research Area

## Image Data



Immunohistochemistry analysis of paraffin-embedded human pancreas, using PEX11C Antibody. The picture on the right is blocked with the synthesized peptide.

## Note

For research use only.