

Product Name: Pxmp2 Rabbit Polyclonal Antibody
Catalog #: APRab16709



Summary

Production Name	Pxmp2 Rabbit Polyclonal Antibody
Description	Rabbit Polyclonal Antibody
Host	Rabbit
Application	IHC,ELISA
Reactivity	Human,Rat,Mouse

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Polyclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% New type preservative N.
Purification	Affinity purification

Immunogen

Gene Name	PXMP2
Alternative Names	PXMP2; PMP22; Peroxisomal membrane protein 2; 22 kDa peroxisomal membrane protein
Gene ID	5827.0
SwissProt ID	Q9NR77.The antiserum was produced against synthesized peptide derived from human PXMP2. AA range:38-87

Application

Dilution Ratio	IHC 1:100-1:300 ELISA: 1:20000
Molecular Weight	

Background


Product Name: Pxmp2 Rabbit Polyclonal Antibody
Catalog #: APRab16709



function: Seems to be involved in pore-forming activity and may contribute to the unspecific permeability of the peroxisomal membrane., similarity: Belongs to the peroxisomal membrane protein PXMP2/4 family., subunit: Interacts with PEX19 and SIVA1., function: Seems to be involved in pore-forming activity and may contribute to the unspecific permeability of the peroxisomal membrane., similarity: Belongs to the peroxisomal membrane protein PXMP2/4 family., subunit: Interacts with PEX19 and SIVA1.,

Research Area

Image Data



Immunohistochemistry analysis of paraffin-embedded human breast carcinoma, using PXMP2 Antibody. The picture on the right is blocked with the synthesized peptide.

Note

For research use only.