

Product Name: Rab L5 Rabbit Polyclonal Antibody
Catalog #: APRab16765



Summary

Production Name	Rab L5 Rabbit Polyclonal Antibody
Description	Rabbit Polyclonal Antibody
Host	Rabbit
Application	IHC,ELISA
Reactivity	Human,Rat,Mouse

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Polyclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% New type preservative N.
Purification	Affinity purification

Immunogen

Gene Name	RABL5
Alternative Names	RABL5; Rab-like protein 5
Gene ID	64792.0
SwissProt ID	Q9H7X7.The antiserum was produced against synthesized peptide derived from human RABL5. AA range:128-177

Application

Dilution Ratio	IHC 1:100-1:300 ELISA: 1:5000
Molecular Weight	

Background

similarity:Belongs to the small GTPase superfamily. Rab family.,similarity:Belongs to the small GTPase superfamily. Rab

