

Product Name: RSAD1 Rabbit Polyclonal Antibody
Catalog #: APRab17400



Summary

Production Name	RSAD1 Rabbit Polyclonal Antibody
Description	Rabbit Polyclonal Antibody
Host	Rabbit
Application	IHC, WB, ELISA
Reactivity	Human, Mouse

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Polyclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% New type preservative N.
Purification	Affinity purification

Immunogen

Gene Name	RSAD1
Alternative Names	RSAD1; Radical S-adenosyl methionine domain-containing protein 1; mitochondrial; Oxygen-independent coproporphyrinogen-III oxidase-like protein RSAD1
Gene ID	55316.0
SwissProt ID	Q9HA92. The antiserum was produced against synthesized peptide derived from human RSAD1. AA range: 298-347

Application

Dilution Ratio	WB 1:500 - 1:2000. IHC 1:100 - 1:300. ELISA: 1:5000..
Molecular Weight	52kD

Product Name: RSAD1 Rabbit Polyclonal Antibody
Catalog #: APRab17400




Background


cofactor: Binds 1 4Fe-4S cluster. The cluster is coordinated with 3 cysteines and an exchangeable S-adenosyl-L-methionine. function: May be involved in porphyrin cofactor biosynthesis. similarity: Belongs to the anaerobic coproporphyrinogen-III oxidase family. cofactor: Binds 1 4Fe-4S cluster. The cluster is coordinated with 3 cysteines and an exchangeable S-adenosyl-L-methionine. function: May be involved in porphyrin cofactor biosynthesis. similarity: Belongs to the anaerobic coproporphyrinogen-III oxidase family.

Research Area

Image Data




Immunohistochemistry analysis of paraffin-embedded human breast carcinoma tissue, using RSAD1 Antibody. The picture on the right is blocked with the synthesized peptide.



Western blot analysis of lysates from HeLa cells, using RSAD1 Antibody. The lane on the right is blocked with the synthesized peptide.

Product Name: RSAD1 Rabbit Polyclonal Antibody
Catalog #: APRab17400



Western Blot analysis of various cells using RSAD1 Polyclonal Antibody

Note

For research use only.