

**Product Name: STRN4 Rabbit Polyclonal Antibody**  
**Catalog #: APRab18410**



## Summary

<b>Production Name</b>	STRN4 Rabbit Polyclonal Antibody
<b>Description</b>	Rabbit Polyclonal Antibody
<b>Host</b>	Rabbit
<b>Application</b>	WB
<b>Reactivity</b>	Human,Mouse

## Performance

<b>Conjugation</b>	Unconjugated
<b>Modification</b>	Unmodified
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>Form</b>	Liquid
<b>Storage</b>	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
<b>Buffer</b>	Liquid in PBS containing 50% glycerol, and 0.02% New type preservative N.
<b>Purification</b>	Affinity purification

## Immunogen

<b>Gene Name</b>	STRN4 ZIN
<b>Alternative Names</b>	
<b>Gene ID</b>	29888.0
<b>SwissProt ID</b>	Q9NRL3.Synthesized peptide derived from human protein . at AA range: 250-330

## Application

<b>Dilution Ratio</b>	WB 1:500-2000 ELISA 1:5000-20000
<b>Molecular Weight</b>	82kD

## Background

function: Binds calmodulin in a calcium dependent manner. May function as scaffolding or signaling

protein.,miscellaneous: The name 'Zinedin' probably originates from the name of the famous soccer player from Marseille


**Product Name: STRN4 Rabbit Polyclonal Antibody**  
**Catalog #: APRab18410**



(Zinedine Zidane),,similarity:Belongs to the WD repeat striatin family,,similarity:Contains 7 WD repeats,,function:Binds calmodulin in a calcium dependent manner. May function as scaffolding or signaling protein,,miscellaneous:The name 'Zinedin' probably originates from the name of the famous soccer player from Marseille (Zinedine Zidane),,similarity:Belongs to the WD repeat striatin family,,similarity:Contains 7 WD repeats,,

## Research Area

## Image Data



Western Blot analysis of HeLa lysis, using primary antibody at 1:1000 dilution. Secondary antibody was diluted at 1:10000

## Note

For research use only.