

Product Name: Anti-Human CD3 (Brenetafusp Biosimilar)(RUO)

Catalog #: ABS00099

For research use only.

Summary

Description	Biosimilar Reference Antibody
Cas.No	2736407-54-6
Source	na
Application	Elisa,FC,Functional assays in vivo
Reactivity	Human
Conjugation	Unconjugated
Isotype	scFv - kappa - heavy - TR-BETA-2_TR-ALPHA
Immunogen	Human CD3
Purity	>98% by SDS-PAGE
Endotoxin level	<1.0 EU/mg as determined by LAL method
Form	Liquid
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	0.01M PBS, pH 7.4