

Product Name: Anti-Human TNFSF11/CD254/RANKL/TRANSE/OPGL/ODF (Denosumab Biosimilar)(RUO)


Catalog #: ABS00177

For research use only.


Summary

Description	Biosimilar Reference Antibody
Cas.No	615258-40-7
Source	Genetically human
Application	Elisa,FC,Functional assays in vivo
Reactivity	Human
Conjugation	Unconjugated
Isotype	IgG2
Immunogen	Human TNFSF11/CD254/RANKL/TRANSE/OPGL/ODF
Purity	>98% by SDS-PAGE
Endotoxin level	<1.0 EU/mg as determined by LAL method
Form	Liquid
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	0.01M PBS, pH 7.4

Image Data



Anti-Human TNFSF11/CD254/RANKL/TRANSE/OPGL/ODF (Denosumab Biosimilar) (ABS00177) on SDS-PAGE under non-reducing (NR) and reducing (R) condition.



ELISA plate pre-coated Antigen protein, can bind Denosumab biosimilar in a linear range of 3.91-125 ng/mL.