


**Product Name:** Anti-Human TIGIT/WUCAM/VSTM3 (Domvanalimab Biosimilar)(RUO)  
**Catalog #:** ABS00189

For research use only.


## Summary

|                        |                                                                             |
|------------------------|-----------------------------------------------------------------------------|
| <b>Description</b>     | Biosimilar Reference Antibody                                               |
| <b>Cas.No</b>          | 2368219-35-4                                                                |
| <b>Source</b>          | Humanised                                                                   |
| <b>Application</b>     | Elisa,FC,Functional assays in vivo                                          |
| <b>Reactivity</b>      | Human                                                                       |
| <b>Conjugation</b>     | Unconjugated                                                                |
| <b>Isotype</b>         | IgG1 - kappa                                                                |
| <b>Immunogen</b>       | Human TIGIT/WUCAM/VSTM3                                                     |
| <b>Purity</b>          | >98% by SDS-PAGE                                                            |
| <b>Endotoxin level</b> | <1.0 EU/mg as determined by LAL method                                      |
| <b>Form</b>            | Liquid                                                                      |
| <b>Storage</b>         | Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles. |
| <b>Shipping</b>        | Ice bags                                                                    |
| <b>Buffer</b>          | 0.01M PBS, pH 7.4                                                           |

## Image Data



Anti-Human TIGIT/WUCAM/VSTM3 (Domvanalimab Biosimilar)(ABS00189) on SDS-PAGE under non-reducing (NR) and reducing (R) condition.



ELISA plate pre-coated Antigen protein, can bind Domvanalimab biosimilar in a linear range of 3.91-125 ng/mL.