

Product Name: Anti-Human TNFSF15/TL1A/VEGI (Duvakitug Biosimilar)(RUO)


Catalog #: ABS00200

For research use only.


Summary

Description	Biosimilar Reference Antibody
Cas.No	2750005-84-4
Source	Genetically human
Application	Elisa,FC,Functional assays in vivo
Reactivity	Human
Conjugation	Unconjugated
Isotype	IgG1 - lambda
Immunogen	Human TNFSF15/TL1A/VEGI
Purity	>98% by SDS-PAGE
Endotoxin level	<1.0 EU/mg as determined by LAL method
Form	Liquid
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	0.01M PBS, pH 7.4

Image Data



Anti-Human TNFSF15/TL1A/VEGI (Duvakitug Biosimilar)(ABS00200) on SDS-PAGE under non-reducing (NR) and reducing (R) condition.



ELISA plate pre-coated Antigen protein, can bind Duvakitug biosimilar in a linear range of 15.63-500 ng/mL.