
Product Name: Anti-Human MS4A1/CD20 (Eramkafusp Biosimilar)(RUO)

Catalog #: ABS00238

For research use only.

Summary

Description	Biosimilar Reference Antibody
Cas.No	2388499-82-7
Source	Chimeric (Murine/Human)
Application	Elisa,FC,Functional assays in vivo
Reactivity	Human
Conjugation	Unconjugated
Isotype	IgG1 - kappa
Immunogen	Human MS4A1/CD20
Purity	>98% by SDS-PAGE
Endotoxin level	<1.0 EU/mg as determined by LAL method
Form	Liquid
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	0.01M PBS, pH 7.4