

**Product Name: Anti-Human IAP/CD47 (Lemzoparlimab Biosimilar)(RUO)**


**Catalog #: ABS00370**

For research use only.


**Summary**

|                        |                                                                             |
|------------------------|-----------------------------------------------------------------------------|
| <b>Description</b>     | Biosimilar Reference Antibody                                               |
| <b>Cas.No</b>          | 2377483-71-9                                                                |
| <b>Source</b>          | Genetically human                                                           |
| <b>Application</b>     | Elisa,FC,Functional assays in vivo                                          |
| <b>Reactivity</b>      | Human                                                                       |
| <b>Conjugation</b>     | Unconjugated                                                                |
| <b>Isotype</b>         | IgG4 - kappa                                                                |
| <b>Immunogen</b>       | Human IAP/CD47                                                              |
| <b>Purity</b>          | >98% by SDS-PAGE                                                            |
| <b>Endotoxin level</b> | <1.0 EU/mg as determined by LAL method                                      |
| <b>Form</b>            | Liquid                                                                      |
| <b>Storage</b>         | Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles. |
| <b>Shipping</b>        | Ice bags                                                                    |
| <b>Buffer</b>          | 0.01M PBS, pH 7.4                                                           |

**Image Data**



Anti-Human IAP/CD47 (Lemzoparlimab Biosimilar)(ABS00370) on SDS-PAGE under non-reducing (NR) and reducing (R) condition.



ELISA plate pre-coated Antigen protein, can bind Lemzoparlimab biosimilar in a linear range of 15.63–125 ng/mL.