

Product Name: Anti-Human TfR1/CD71 (Lepunafusp Biosimilar)(RUO)

Catalog #: ABS00373

For research use only.

Summary

Description	Biosimilar Reference Antibody
Cas.No	2396411-71-3
Source	na
Application	Elisa,FC,Functional assays in vivo
Reactivity	Human
Conjugation	Unconjugated
Isotype	F(ab') ₂ - G1 - kappa
Immunogen	Human TfR1/CD71
Purity	>98% by SDS-PAGE
Endotoxin level	<1.0 EU/mg as determined by LAL method
Form	Liquid
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	0.01M PBS, pH 7.4