

Product Name: Anti-Human LY75/CD205/DEC205 (Oberotatug Biosimilar)(RUO)


Catalog #: ABS00482

For research use only.


Summary

Description	Biosimilar Reference Antibody
Cas.No	2641919-35-7
Source	Genetically human
Application	Elisa,FC,Functional assays in vivo
Reactivity	Human
Conjugation	Unconjugated
Isotype	IgG1 - kappa
Immunogen	Human LY75/CD205/DEC205
Purity	>98% by SDS-PAGE
Endotoxin level	<1.0 EU/mg as determined by LAL method
Form	Liquid
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	0.01M PBS, pH 7.4

Image Data



Anti-Human LY75/CD205/DEC205 (Oberotatug Biosimilar)(ABS00482) on SDS-PAGE under non-reducing (NR) and reducing (R) condition.



ELISA plate pre-coated Antigen protein, can bind Oberotatug biosimilar in a linear range of 3.91–31 ng/mL.