

Product Name: Anti-Human PDCD1/CD279/PD1 (Penpulimab Biosimilar)(RUO)


Catalog #: ABS00539

For research use only.


Summary

Description	Biosimilar Reference Antibody
Cas.No	2350298-92-7
Source	Chimeric and/or humanised
Application	Elisa,FC,Functional assays in vivo
Reactivity	Human
Conjugation	Unconjugated
Isotype	IgG1 - kappa
Immunogen	Human PDCD1/CD279/PD1
Purity	>98% by SDS-PAGE
Endotoxin level	<1.0 EU/mg as determined by LAL method
Form	Liquid
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	0.01M PBS, pH 7.4

Image Data



Anti-Human PDCD1/CD279/PD1 (Penpulimab Biosimilar)(ABS00539) on SDS-PAGE under non-reducing (NR) and reducing (R) condition.



ELISA plate pre-coated Antigen protein, can bind Penpulimab biosimilar in a linear range of 0.18–15 ng/mL.