


Product Name: Anti-Human PDCD1/CD279/PD1 (Retifanlimab Biosimilar)(RUO)**Catalog #: ABS00602**


For research use only.

Summary

Description	Biosimilar Reference Antibody
Cas.No	2079108-44-2
Source	Chimeric and Humanised
Application	Elisa,FC,Functional assays in vivo
Reactivity	Human
Conjugation	Unconjugated
Isotype	IgG4 - kappa
Immunogen	Human PDCD1/CD279/PD1
Purity	>98% by SDS-PAGE
Endotoxin level	<1.0 EU/mg as determined by LAL method
Form	Liquid
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	0.01M PBS, pH 7.4

Image Data

Anti-Human PDCD1/CD279/PD1 (Retifanlimab Biosimilar)(ABS00602) on SDS-PAGE under non-reducing (NR) and reducing (R) condition.



ELISA plate pre-coated Antigen protein, can bind Retifanlimab biosimilar in a linear range of 15.63–125 ng/mL.