
Product Name: Anti-Human IL36RN (Spesolimab Biosimilar)(RUO)

Catalog #: ABS00675

For research use only.

Summary

Description	Biosimilar Reference Antibody
Cas.No	2097104-58-8
Source	Humanised
Application	Elisa,FC,Functional assays in vivo
Reactivity	Human
Conjugation	Unconjugated
Isotype	IgG1 - kappa
Immunogen	Human IL36RN
Purity	>98% by SDS-PAGE
Endotoxin level	<1.0 EU/mg as determined by LAL method
Form	Liquid
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	0.01M PBS, pH 7.4