

Product Name: Anti-Canine IL31 (Canine) (Tirnovetmab Biosimilar)(RUO)


Catalog #: ABS00728

For research use only.

Summary

Description	Biosimilar Reference Antibody
Cas.No	2364504-80-1
Source	Caninised
Application	Elisa,FC,Functional assays in vivo
Reactivity	Canine
Conjugation	Unconjugated
Isotype	IgG2 - kappa
Immunogen	Canine IL31
Purity	>98% by SDS-PAGE
Endotoxin level	<1.0 EU/mg as determined by LAL method
Form	Liquid
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	0.01M PBS, pH 7.4

Image Data



Anti-Canine IL31 (Canine) (Tirnovetmab Biosimilar)(ABS00728) on SDS-PAGE under non-reducing (NR) and reducing (R) condition.