

Ham's F-10 Complete Medium Product manual

Basic Information

Cat.NO	Size	Shelf	Form	Storage	Transportation
CMC0045	500mL	6 months	Liquid	Store at 2-8°C away from light	2-8°C

Product Introduction

Components	Ham's F-10(CMB0086)+10% Nutrient+1% Supplement
Concentration	Ready-to-use
Bacteria Detection	Negative
Fungi Detection	Negative
Mycoplasma Detection	Negative
Endotoxin Content (EU/mL)	< 3

Background

Ham's F-10 Nutrient Mixture is a culture medium designed by Ham in 1963 for serum-free culture of CHO cells. Ham's F-10 contains rich nutrients and is suitable for the culture of human diploid cells, cell cloning culture under low serum conditions, normal culture of mammalian cells, chromosome analysis of leukocytes, and primary cell culture of rats, rabbits, chickens and other animals.

Precautions

1. During the entire process, be sure to pay attention to aseptic operation to avoid contamination;
2. To maintain the best use effect of this product, please do not perform freeze-thaw treatment;
3. This product is only used for research or further research, not for diagnosis and treatment.
4. It is not suitable for long time storage at room temperature.
5. This product is a ready-to-use medium. If there is no special need, don't add serum, penicillin and streptomycin. It can be used directly.