

## Product Information

Our pathology-grade Antibody Diluent is a purpose-formulated buffer for diluting diagnostic primary antibodies. Optimized for both IHC and IF applications, it incorporates protective stabilizers that preserve antibody activity and extend shelf life while suppressing non-specific binding, cutting background staining, and boosting assay specificity and sensitivity. Choose our pathology-dedicated diluent for sharper, faster, and more reliable tissue-based diagnostics—empowering your research with every slide.

## Basic Information

Product Form	liquid
Size	50mL/100mL/500mL
Storage	2-8°C
Transportation	Packed with Ice bag, 2-8°C
Validity	12 months

## Operation steps

1. Dilute the antibody concentrate to the recommended working concentration with Antibody Diluent.
2. Directly add antibody working solution onto tissue sections or experimental samples.

## Notes

1. For your safety and health, please wear a lab coat, disposable gloves and mask during operation.
2. Improper storage of the antibody diluent may lead to protein degradation and loss of activity, always store under the conditions specified in the instructions.
3. This product is For Research Use Only, Not for Diagnostic Use.