
Product Name: Prolactin(PRL) Mouse Monoclonal Antibody**Catalog #: AMM21980**

For research use only.

Summary

Description	Mouse Monoclonal Antibody
Host	Mouse
Application	IHC,IF,ELISA
Reactivity	Human
Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG1,Kappa
Clonality	Monoclonal
Form	Liquid
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	PBS, 50% glycerol, 0.05% Proclin 300, 0.05%BSA
Purification	The antibody was affinity-purified from ascites by affinity-chromatography using specific immunogen.

Application

Dilution Ratio	IHC 1:200-400;IF 1:50-200;ELISA 1:500-5000
Molecular Weight	Calculated MW:26kDa

Antigen Information

Gene Name	PRL
Alternative Names	
Gene ID	Human:5617
SwissProt ID	Human:P01236
Immunogen	Synthesized peptide derived from human Prolactin AA range: 100-227

Background


This gene encodes the anterior pituitary hormone prolactin. This secreted hormone is a growth regulator for many tissues, including cells of the immune system. It may also play a role in cell survival by suppressing apoptosis, and it is essential for lactation. Alternative splicing results in multiple transcript variants that encode the same protein. [provided by RefSeq, Aug

2011],

Research Area

Pathology

Image Data



Human pituitary adenoma tissue was stained with Anti-Prolactin Antibody