
Product Name: CA IX Mouse Monoclonal Antibody**Catalog #: AMM22043**

For research use only.

Summary

Description	Mouse Monoclonal Antibody
Host	Mouse
Application	IHC,ELISA
Reactivity	Human
Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	PBS, 50% glycerol, 0.05% Proclin 300, 0.05%BSA
Purification	The antibody was affinity-purified from ascites by affinity-chromatography using specific immunogen.

Application

Dilution Ratio	IHC 1:200-400;ELISA 1:500-5000
Molecular Weight	Calculated MW:49kDa

Antigen Information


Gene Name	CA9 G250;MN;Carbonic anhydrase 9;Carbonate dehydratase IX;Carbonic anhydrase IX;CA-
Alternative Names	IX;CAIX;Membrane antigen MN;P54/58N;Renal cell carcinoma-associated antigen G250;RCC-associated antigen G250;pMW1
Gene ID	Human:768,Mouse:230099
SwissProt ID	Human:Q16790,Mouse:Q8VHB5
Immunogen	Recombinant protein

Background

Research Area

Pathology

Image Data



Human stomach tissue was stained with Anti-CA IX Antibody