
Product Name: Cytokeratin 19(CK19) Mouse Monoclonal Antibody**Catalog #: AMM22045**

For research use only.

Summary

Description	Mouse Monoclonal Antibody
Host	Mouse
Application	IHC,ELISA
Reactivity	Human
Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG1,Kappa
Clonality	Monoclonal
Form	Liquid
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	PBS, 50% glycerol, 0.05% Proclin 300, 0.05%BSA
Purification	Protein G

Application

Dilution Ratio	IHC 1:200-400.ELISA 1:500-5000.
Molecular Weight	Calculated MW:44kDa,Observed MW:44kDa

Antigen Information

Gene Name	KRT19
Alternative Names	Keratin, type I cytoskeletal 19;Cytokeratin-19;CK-19;Keratin-19;K19;
Gene ID	Human:3880
SwissProt ID	Human:P08727
Immunogen	Recombinant protein

Background


The protein encoded by this gene is a member of the keratin family. The keratins are intermediate filament proteins responsible for the structural integrity of epithelial cells and are subdivided into cytokeratins and hair keratins. The type I cytokeratins consist of acidic proteins which are arranged in pairs of heterotypic keratin chains. Unlike its related family members, this smallest known acidic cytokeratin is not paired with a basic cytokeratin in epithelial cells. It is specifically expressed in the

periderm, the transiently superficial layer that envelops the developing epidermis. The type I cytokeratins are clustered in a region of chromosome 17q12-q21. [provided by RefSeq, Jul 2008],

Research Area

Pathology

Image Data



Human appendix tissue was stained with Anti-Cytokeratin 19 Antibody