
Product Name: Melan-A(Mart1) mouse Monoclonal Antibody**Catalog #: AMM22083**

For research use only.

Summary

Description	Mouse Monoclonal Antibody
Host	Mouse
Application	IHC,ELISA
Reactivity	Human
Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG2b,Kappa
Clonality	Monoclonal
Form	Liquid
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	PBS, 50% glycerol, 0.05% Proclin 300, 0.05%BSA
Purification	The antibody was affinity-purified from ascites by affinity-chromatography using specific immunogen.

Application

Dilution Ratio	IHC 1:200-400;ELISA 1:500-5000
Molecular Weight	Calculated MW:13kDa,Observed MW:18kDa

Antigen Information

Gene Name	MLANA MART1
Alternative Names	Melanoma antigen recognized by T-cells 1;MART-1;Antigen LB39-AA;Antigen SK29-AA;Protein Melan-A;
Gene ID	Human:2315
SwissProt ID	Human:Q16655
Immunogen	Synthesized peptide derived from human Mart-1/Melan-A AA range: 50-118


Background

tissue specificity:Expression is restricted to melanoma and melanocyte cell lines and retina.,

Research Area

Pathology

Image Data



Human malignant mesothelioma tissue was stained with Anti-Melan-A Antibody