
Product Name: Cadherin Mouse Monoclonal Antibody**Catalog #: AMM22099**

For research use only.

Summary

Description	Mouse Monoclonal Antibody
Host	Mouse
Application	IHC,ELISA
Reactivity	Human,Mouse,Rat
Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG2a,Kappa
Clonality	Monoclonal
Form	Liquid
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	PBS, 50% glycerol, 0.05% Proclin 300, 0.05%BSA
Purification	The antibody was affinity-purified from ascites by affinity-chromatography using specific immunogen.

Application

Dilution Ratio	IHC 1:200-400;ELISA 1:500-5000
Molecular Weight	Observed MW:130kDa

Antigen Information


Gene Name	cadherin (ABT123)
Alternative Names	
Gene ID	
SwissProt ID	Human:P55285;P22223;P12830;P55283
Immunogen	Synthesized peptide derived from human cadherin AA range: 800-882 according to Cadherin-1

Background

Research Area

Pathology

Image Data



Human appendix tissue was stained with anti-cadherin antibody.