
Product Name: Mammaglobin Mouse Monoclonal Antibody**Catalog #: AMM22132**

For research use only.

Summary

Description	Mouse Monoclonal Antibody
Host	Mouse
Application	IHC,IF,ELISA
Reactivity	Human
Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG1,Kappa
Clonality	Monoclonal
Form	Liquid
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	PBS, 50% glycerol, 0.05% Proclin 300, 0.05%BSA
Purification	The antibody was affinity-purified from ascites by affinity-chromatography using specific immunogen.

Application

Dilution Ratio IHC 1:200-400;IF 1:50-200;ELISA 1:500-5000

Molecular Weight


Antigen Information

Gene Name	SCGB2A2 MGB1 UGB2
Alternative Names	Mammaglobin-A;Mammaglobin-1;Secretoglobin family 2A member 2;
Gene ID	Human:4250
SwissProt ID	Human:Q13296;O75556
Immunogen	Synthesized peptide derived from human Mammaglobin AA range: 19-93

Background**Research Area**

Pathology

Image Data



Human breast tissue was stained with Anti-Mammaglobin Antibody