



Product Name: Recombinant Human CS Protein
Catalog #: PBH2743



Summary

Name	CISY_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	O75390
Host	Baculovirus infected Sf9
Species	Human
Molecular Mass	53 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	MALLTAAARL LGTKNASCLV LAARHASASS TNLKDILADL IPKEQARIKT FRQQHGKTVV GQITVDMMYG GMRGMKGLVY ETSVLDPDEG IRFRGFSIPE CQKLLPKAKG GEEPLPEGLF WLLVTGHIPT EEQVSWLSKE WAKRAALPSH VVTMLDNFPT NLHPMSQLSA AVTALNSESF FARAYAQQIS RTKYWELIYE DSMDLIAKLP CVAAKIYRNL YREGSGIGAI DSNLDWSHNF TNMLGYTDHQ FTELTRYLYT IHSDEHGGNV SAHTSHLVGS ALSDPYLSFA AAMNGLAGPL HGLANQEVLV WLTQLQKEVG KDVSEKLRD YIWNTLNLSGR VVPGYGHAVL RKTDPRYTCQ REFALKHLPN DPMFKLVAQL YKIVPNVLE QGKAKNPWPN VDAHSGVLLQ YGGMTEMNYY TVLFGVSRAL GVLAQLIWSR ALGFPLERPK SMSTEGLMKF VDSKSG

image



Background

Product Name: Recombinant Human CS Protein
Catalog #: PBH2743



Alternative Names

Citrate synthase, mitochondrial (EC 2.3.3.1) (Citrate (Si)-synthase)

Note

For Research Use Only , Not for Diagnostic Use