



Product Name: Recombinant Human TPP1 Protein
Catalog #: PBH5563



Summary

Name	TPP1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	O14773
Host	Baculovirus infected Sf9
Species	Human
Molecular Mass	62 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	MGLQACLLGL FALISGKCS YSPEPDQRRT LPPGWWSLGR ADPEEELSLT FALRQQNVER LSELVQAVSD PSSPQYGKYL TLENVADLVR PSPLTLHTVQ KWLLAAGAQK CHSVITQDFL TCWLSIRQAE LLLPGAEFHH YVGGPTETHV VRSPHPYQLP QALAPHVDFV GGLHRFPPTS SLRQRPEPQV TGTVGLHLGV TPSVIRKRYN LTSQDVSGT SNNSQACAQF LEQYFHSDSL AQFMRLFGGN FAHQASVARV VGQQGRGRAG IEASLDVQYL MSAGANISTW VYSSPGRHEG QEPFLQWLML LSNESALPHV HTVSYGDDSD SLSSAYIQRV NTELMKAAAR GLTLLFASGD SGAGCWSVSG RHQFRPTFPA SSPYVTTVGG TSFQEPFLIT NEIVDYISGG GFSNVFPRPS YQEEAVTKFL SSSPHLPPSS YFNASGRAYP DVAALSDGYW VVSNRVPIPW VSGTSASTPV FGGILSLINE HRILSGRPPL GFLNPRLYQQ HGAGLFDVTR GCHESCLDEE VEGQGFCSGP GWDPVTGWGT PNFPALLKTL LNP

image



Product Name: Recombinant Human TPP1 Protein
Catalog #: PBH5563



Background

Alternative Names

Tripeptidyl-peptidase 1 (TPP-1) (EC 3.4.14.9) (Cell growth-inhibiting gene 1 protein) (Lysosomal pepstatin-insensitive protease) (LPIC) (Tripeptidyl aminopeptidase) (Tripeptidyl-peptidase I) (TPP-I)

Note

For Research Use Only , Not for Diagnostic Use