


Product Name: Recombinant Human IgE Protein
Catalog #: PCH10924



Summary

Name	IGHE_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	tag free
Accession #	P01854
Host	CHO
Species	Human
Molecular Mass	kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	序列

image



Background

Alternative Names	Immunoglobulin heavy constant epsilon (Ig epsilon chain C region) (Ig epsilon chain C region ND)
--------------------------	--------------------------------------------------------------------------------------------------

Note

For Research Use Only , Not for Diagnostic Use