

Product Name: Recombinant Human ANXA10
Catalog #: PEH0069




Summary

Name	Annexin A10/ANXA10
Purity	Greater than 95% as determined by reducing SDS-PAGE
Endotoxin level	<1 EU/μg as determined by LAL test.
Construction	Recombinant Human Annexin A10 is produced by our E.coli expression system and the target gene encoding Met1-Tyr324 is expressed.
Accession #	Q9UJ72
Host	E.coli
Species	Human
Predicted Molecular Mass	37.5 KDa
Formulation	Lyophilized from a 0.2 μm filtered solution of PBS, pH7.4.
Shipping	The product is shipped at ambient temperature. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Lyophilized protein should be stored at ≤ -20°C, stable for one year after receipt. Reconstituted protein solution can be stored at 2-8°C for 2-7 days. Aliquots of reconstituted samples are stable at ≤ -20°C for 3 months.
Reconstitution	Always centrifuge tubes before opening. Do not mix by vortex or pipetting. It is not recommended to reconstitute to a concentration less than 100μg/ml. Dissolve the lyophilized protein in distilled water. Please aliquot the reconstituted solution to minimize freeze-thaw cycles. Always centrifuge tubes before opening. Do not mix by vortex or pipetting. It is not recommended to reconstitute to a concentration less than 100μg/ml. Dissolve the lyophilized protein in distilled water. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.

SDS-PAGE image

Product Name: Recombinant Human ANXA10
Catalog #: PEH0069



Alternative Names

Annexin A10; Annexin-10; Annexin-14; ANXA10; ANX14

Background

Annexin A10 (ANXA10) contains four Annexin repeats and is a member of the Annexin family. Members of this calcium-dependent phospholipid-binding protein family play a role in the regulation of cellular growth and in signal transduction pathways. It is reported that ANXA10 may be a clinical relevant marker for predicting outcome in both early and advanced stages of bladder cancer.

Note

For Research Use Only , Not for Diagnostic Use.