

**Product Name: Recombinant Human NT-proBNP (N-6His-Flag)**  
**Catalog #: PEH0167**




---

## Summary

<b>Name</b>	NT-proBNP
<b>Purity</b>	Greater than 95% as determined by reducing SDS-PAGE
<b>Endotoxin level</b>	<1 EU/μg as determined by LAL test.
<b>Construction</b>	Recombinant Human Brain-type Natriuretic Peptide is produced by our E.coli expression system and the target gene encoding His27-Arg102 is expressed with a 6His, Flag tag at the N-terminus.
<b>Accession #</b>	P16860
<b>Host</b>	E.coli
<b>Species</b>	Human
<b>Predicted Molecular Mass</b>	11 KDa
<b>Formulation</b>	Lyophilized from a 0.2 μm filtered solution of 20mM Tris-HCl, 150mM NaCl, pH 8.0.
<b>Shipping</b>	The product is shipped at ambient temperature. Upon receipt, store it immediately at the temperature listed below.
<b>Stability&amp;Storage</b>	Lyophilized protein should be stored at ≤ -20°C, stable for one year after receipt. Reconstituted protein solution can be stored at 2-8°C for 2-7 days. Aliquots of reconstituted samples are stable at ≤ -20°C for 3 months.
<b>Reconstitution</b>	Always centrifuge tubes before opening. Do not mix by vortex or pipetting. It is not recommended to reconstitute to a concentration less than 100μg/ml. Dissolve the lyophilized protein in distilled water. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.

## SDS-PAGE image

**Product Name: Recombinant Human NT-proBNP (N-6His-Flag)**  
**Catalog #: PEH0167**



### **Alternative Names**

Natriuretic peptides B; Gamma-brain natriuretic peptide; NPPB; BNP

### **Background**

Brain-type Natriuretic Peptide (BNP) is a nonglycosylated peptide that is produced predominantly by ventricular myocytes and belongs to the natriuretic peptide family. Proteolytic cleavage of the 12 kDa BNP precursor gives rise to N-terminal Pro BNP (NT-proBNP) and mature BNP. N-terminal proB-type natriuretic peptide (NT-proBNP), a useful marker of heart failure (HF), is considered to be secreted mainly from the ventricle, increased serum NT-proBNP levels are also encountered in conditions such as atrial fibrillation (AF) and atrial septal defect in patients without HF.

### **Note**

For Research Use Only , Not for Diagnostic Use.