


Product Name: Recombinant Human CRYAB (C-6His)
Catalog #: PEH0454



Summary

Name	CRYAB/Alpha-crystallin B chain
Purity	Greater than 95% as determined by reducing SDS-PAGE
Endotoxin level	<1 EU/μg as determined by LAL test.
Construction	Recombinant Human Crystalline Alpha -B Chain is produced by our E.coli expression system and the target gene encoding Met1-Lys175 is expressed with a 6His tag at the C-terminus.
Accession #	P02511
Host	E.coli
Species	Human
Predicted Molecular Mass	21.22 KDa
Formulation	Lyophilized from a 0.2 μm filtered solution of PBS, pH 7.4.
Shipping	The product is shipped at ambient temperature. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at ≤-70°C, stable for 6 months after receipt. Store at ≤-70°C, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Reconstitution	Always centrifuge tubes before opening. Do not mix by vortex or pipetting. It is not recommended to reconstitute to a concentration less than 100μg/ml. Dissolve the lyophilized protein in distilled water. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.

SDS-PAGE image



Background

Product Name: Recombinant Human CRYAB (C-6His)
Catalog #: PEH0454



Alternative Names

Alpha-Crystallin B Chain; Alpha(B)-Crystallin; Heat Shock Protein Beta-5; HspB5; Renal Carcinoma Antigen NY-REN-27; Rosenthal Fiber Component; CRYAB; CRYA2

Background

α Crystallin B Chain (CRYAB) is a cytoplasmic protein that belongs to the small heat shock protein (HSP20) family. Alpha crystallins are composed of two gene products: alpha-A and alpha-B, for acidic and basic, respectively. Alpha crystallins can be induced by heat shock and are members of the small heat shock protein (sHSP also known as the HSP20) family. Alpha crystallins acts as molecular chaperones and hold them in in large soluble aggregates. CRYAB is expressed widely in many tissues and organs. It may contribute to the transparency and refractive index of the lens. The deficiency of CRYAB is the cause of myopathy myofibrillar type 2 (MFM2) and cataract posterior polar type 2 (CTPP2).

Note

For Research Use Only , Not for Diagnostic Use.