

Product Name: Recombinant Human FIS1 (C-6His)
Catalog #: PEH0668




Summary

Name	FIS1/Mitochondrial fission 1 protein/TTC11/CGI-135
Purity	Greater than 95% as determined by reducing SDS-PAGE
Endotoxin level	<1 EU/μg as determined by LAL test.
Construction	Recombinant Human Mitochondrial Fission 1 Protein is produced by our E.coli expression system and the target gene encoding Met1-Gly122 is expressed with a 6His tag at the C-terminus.
Accession #	Q9Y3D6
Host	E.coli
Species	Human
Predicted Molecular Mass	15.3 KDa
Formulation	Lyophilized from a 0.2 μm filtered solution of 20mM Tris-HCl, pH 8.0.
Shipping	The product is shipped at ambient temperature. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Lyophilized protein should be stored at ≤ -20°C, stable for one year after receipt. Reconstituted protein solution can be stored at 2-8°C for 2-7 days. Aliquots of reconstituted samples are stable at ≤ -20°C for 3 months.
Reconstitution	Always centrifuge tubes before opening. Do not mix by vortex or pipetting. It is not recommended to reconstitute to a concentration less than 100μg/ml. Dissolve the lyophilized protein in distilled water. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.

SDS-PAGE image

Product Name: Recombinant Human FIS1 (C-6His)
Catalog #: PEH0668



Alternative Names

Mitochondrial Fission 1 Protein; FIS1 Homolog; hFis1; Tetratricopeptide Repeat Protein 11; TPR Repeat Protein 11; FIS1; TTC11; CGI-135

Background

Mitochondrial Fission 1 Protein (FIS1) is a member of the FIS1 family. FIS1 is a single-pass membrane protein and contains one TPR repeat. FIS1 is part of the mitochondrial complex that promotes mitochondrial fission. FIS1 can induce cytochrome C discharge from the mitochondrion to the cytosol, eventually leading to apoptosis. In addition, FIS1 participates in peroxisomal growth and division. The C-terminus of FIS1 is required for mitochondrial or peroxisomal localization, while the N-terminus is necessary for mitochondrial or peroxisomal fission, localization and regulation of the interaction with DNM1L.

Note

For Research Use Only , Not for Diagnostic Use.