

Product Name: Recombinant Human GADD45G (N-6His)
Catalog #: PEH0699




Summary

Name	GADD45G/CR6/DDIT2
Purity	Greater than 95% as determined by reducing SDS-PAGE
Endotoxin level	<1 EU/μg as determined by LAL test.
Construction	Recombinant Human Growth Arrest And DNA Damage-Inducible Protein GADD45 Gamma is produced by our E.coli expression system and the target gene encoding Met1-Glu159 is expressed with a 6His tag at the N-terminus.
Accession #	O95257
Host	E.coli
Species	Human
Predicted Molecular Mass	19.28 KDa
Formulation	Lyophilized from a 0.2 μm filtered solution of 20mM Tris-HCl, 150mM NaCl, pH 8.0.
Shipping	The product is shipped at ambient temperature. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at ≤-70°C, stable for 6 months after receipt. Store at ≤-70°C, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Reconstitution	Always centrifuge tubes before opening. Do not mix by vortex or pipetting. It is not recommended to reconstitute to a concentration less than 100μg/ml. Dissolve the lyophilized protein in distilled water. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.

SDS-PAGE image

Product Name: Recombinant Human GADD45G (N-6His)
Catalog #: PEH0699



Alternative Names

Growth Arrest and DNA Damage-Inducible Protein GADD45 Gamma; Cytokine-Responsive Protein CR6; DNA Damage-Inducible Transcript 2 Protein; DDIT-2; GADD45G; CR6; DDIT2

Background

Growth Arrest and DNA Damage-Inducible Protein GADD45 Y (GADD45G) is a nuclear protein which belongs to the GADD45 family. GADD45G is highly expressed in placenta. GADD45G interacts with various proteins whose transcript levels are increased following stressful growth arrest conditions and treatment with DNA-damaging agents. GADD45G responds to environmental stresses by mediating activation of the p38/JNK pathway via MTK1/MEKK4 kinase. GADD45G is also involved in the regulation of growth and apoptosis. GADD45G inhibits cell growth and differentiation by androgens. The mRNA expression is down-regulated in hepatocellular carcinoma.

Note

For Research Use Only , Not for Diagnostic Use.