


Product Name: Recombinant Human HSF2 (N-6His)
Catalog #: PEH0808



Summary

Name	HSF2
Purity	Greater than 95% as determined by reducing SDS-PAGE
Endotoxin level	<1 EU/μg as determined by LAL test.
Construction	Recombinant Human Heat Shock Factor Protein-2 is produced by our E.coli expression system and the target gene encoding Ser411-Ser536 is expressed with a 6His tag at the N-terminus.
Accession #	Q03933
Host	E.coli
Species	Human
Predicted Molecular Mass	15.9 KDa
Formulation	Lyophilized from a 0.2 μm filtered solution of 20mM PB, 150mM NaCl, 1mM DTT, pH 7.2.
Shipping	The product is shipped at ambient temperature. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at ≤-70°C, stable for 6 months after receipt. Store at ≤-70°C, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Reconstitution	Always centrifuge tubes before opening. Do not mix by vortex or pipetting. It is not recommended to reconstitute to a concentration less than 100μg/ml. Dissolve the lyophilized protein in distilled water. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.

SDS-PAGE image



Product Name: Recombinant Human HSF2 (N-6His)
Catalog #: PEH0808



Background

Alternative Names

Heat Shock Factor Protein 2; HSF 2; Heat Shock Transcription Factor 2; HSTF 2; HSF2; HSTF2

Background

Heat Shock Factor Protein 2 (HSF2) belongs to the HSF family of transcription factors that bind specifically to the heat-shock promoter element and activate transcription. In higher eukaryotes, HSF is unable to bind to the HSE unless the cells are heat shocked. HSF2 is widely expressed in many cells and tissues. HSF2 is located on Cytoplasmic during normal growth. But when it is activated, HSF2 moves to the nucleus.

Note

For Research Use Only , Not for Diagnostic Use.