


Product Name: Recombinant Human ZBED4 Protein
Catalog #: PEH10000



Summary

Name	ZBED4_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	O75132
Host	E.coli
Species	Human
Molecular Mass	28 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	KFKIEEEDDD GIPDSLERM DFKSEQEDMK QTDSGGERAG LGGTGCSCKP PGKYLSAESE DDYGALFSQY SSTLYDVAME AVTQSLSSR NMSSRKKSPA WKHFFISPRD STKAICMYCV KEFSRGKNEK DLSTSCLMRH VRRAHPTVLI QENGSVSAVS SFPSPSLLLP PQPADAGDLS TILSPIKLVQ KVASKIPSPD RITEESVSVV SSEEISSDMS VSEKCGREEA LVGSSPHLPA

image



Background

Alternative Names Zinc finger BED domain-containing protein 4

Note

Product Name: Recombinant Human ZBED4 Protein
Catalog #: PEH10000



For Research Use Only , Not for Diagnostic Use