


Product Name: Recombinant Human ZBED5 Protein
Catalog #: PEH10001



Summary

Name	ZBED5_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q49AG3
Host	E.coli
Species	Human
Molecular Mass	18 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	NDQKVGILQS EDKQLQPSVS KKSEGELSRV KFISNSNKIT FSKKPKRRKY DESYLSFGFT YFGNRDAPHA QCVLCKKILS NSSLAPSKLR RHLETKHAAY KDKDISFFKQ HLDSPENNKP PTPKIVNTDN ESATEASYNV

image



Background

Alternative Names Zinc finger BED domain-containing protein 5 (Transposon-derived Buster1 transposase-like protein)

Note

For Research Use Only , Not for Diagnostic Use