


Product Name: Recombinant Human ZC2HC1C Protein
Catalog #: PEH10002



Summary

Name	ZC21C_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q53FD0
Host	E.coli
Species	Human
Molecular Mass	15 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	MYKPIFSPEF EEEEEFSRDR REDETWGRSQ QNSGPFQFSD YRIQRLKRER LVASNNKIRD PVSEPSVEKF SPPSETPVGA LQGSARNSSL SMAPDSSGSS GSIEEPQLGE

image



Background

Alternative Names Zinc finger C2HC domain-containing protein 1C

Note

For Research Use Only , Not for Diagnostic Use