


Product Name: Recombinant Human ZC4H2 Protein
Catalog #: PEH10003



Summary

Name	ZC4H2_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9NQZ6
Host	E.coli
Species	Human
Molecular Mass	20 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	LKAEFEALES EERHLKEYKQ EMDLLLQEKM AHVEELRLIH ADINVMENTI KQSENDLNKL LESTRRLHDE YKPLKEHVDA LRMTLGLQRL PDLCEEEEL SLDYFEKQKA EWQTEPQEPP IPESLAAAAA AAQLQVARK QDTRQTATFR

image



Background

Alternative Names Zinc finger C4H2 domain-containing protein (Hepatocellular carcinoma-associated antigen 127)

Note

For Research Use Only , Not for Diagnostic Use