


Product Name: Recombinant Human ZC3H14 Protein
Catalog #: PEH10004



Summary

Name	ZC3HE_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q6PJT7
Host	E.coli
Species	Human
Molecular Mass	26 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	YGSRTGSISS SVSVPAKPER RPSLPPSKQA NKNLILKAIS EAQESVTKTT NYSTVPQKQT LPVAPRTRTS QEELLAEEVVQ GQSRTPRISP PIKEEETKGD SVEKNQGTQQ RQLLSRLQID PVMAETLQMS QDYYDMESMV HADTRSFILK KPKLSEEVVV APNQESGMKT ADSLRVLSGH LMQTRDLVQP DKPASPKFIV TLDGVPSPPG

image



Background

Alternative Names	Zinc finger CCCH domain-containing protein 14 (Mammalian suppressor of tau pathology-2) (MSUT-2) (Renal carcinoma antigen NY-REN-37)
--------------------------	--

Note

Product Name: Recombinant Human ZC3H14 Protein
Catalog #: PEH10004



For Research Use Only , Not for Diagnostic Use