


Product Name: Recombinant Human ZMYND10 Protein
Catalog #: PEH10010



Summary

Name	ZMY10_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	O75800
Host	E.coli
Species	Human
Molecular Mass	33 kDa
Buffer	Lyophilized from Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	EAELVLRGLR SFPLREMGSE GWNQQHENLE KLNMQAILDA TVSQGEPIQE LLVTHGKVPT LVEELIAVEM WKQKVPVFC RVEDFKPQNT FPIYMVVHHE ASIINLLETV FFHKEVCESA EDTVLDLVDY CHRKLTLVA QSGCGGPPEG EGSQDSNPMQ ELQKQAEELME FEALKALSV LRYITDCVDS LSLSTLSRML STHNLPCLLV ELLEHSPWSR REGGKLQQFE GSRWHTVAPS EQQKLSKLDG QVWIALYNLL LSPEAQARYC

image



Background

Alternative Names Zinc finger MYND domain-containing protein 10 (Protein BLu)

Product Name: Recombinant Human ZMYND10 Protein
Catalog #: PEH10010



Note

For Research Use Only , Not for Diagnostic Use